


Wearable Clap Switch

ASSEMBLY GUIDE

Required Components



BOM

- | | |
|---------------------------------|------------------------|
| 1. Transmitter Lid | 18. Female header pins |
| 2. Transmitter Base | 19. M2 Screw x8 |
| 3. Transmitter Cover x2 | 20. Toggle switch |
| 4. Transmitter Button x2 | 21. 5V relay |
| 5. M3 screw x6 | 22. Potentiometer x2 |
| 6. Double Sided Velcro | 23. Breadboard wire |
| 7. Circuit Playground Bluefruit | 24. Wire |
| 8. LiPo Battery | 25. M3 screw x2 |
| 9. Battery Charger | |
| 10. Slide Switch | |
| 11. Receiver Base | |
| 12. Receiver Lid | |
| 13. ItsyBitsy nrf52408 | |
| 14. LiPo Battery | |
| 15. Battery Backpack | |
| 16. 3.5mm mono jack | |
| 17. Protoboard | |

Wearable Clap Switch

ASSEMBLY GUIDE

Required Tools

- Soldering iron
- Flush Cutters
- Wire stripper
- Solder
- Electrical Tape
- Screwdriver

Required Personal Protective Equipment (PPE)

- Safety Glasses

Wearable Clap Switch

ASSEMBLY GUIDE

Assembly Instructions

Transmitter

Step 1


Take two 5cm pieces of wire and strip 1cm from each end.



Wearable Clap Switch ASSEMBLY GUIDE

Step 2


Cut the black wire on the battery 4cm from the JST connector and strip 1cm on each end of the cut.



Wearable Clap Switch ASSEMBLY GUIDE

Step 3

Solder the wires from step 1 onto each end of the stripped cut from step 2 and cover the soldered joints with electrical tape.



Wearable Clap Switch ASSEMBLY GUIDE

Step 4

Solder the remaining two ends of the wire to the middle and side terminal of the slide switch.



Wearable Clap Switch

ASSEMBLY GUIDE

Step 5

Cut two lengths of double sided hook and loop fastener, roughly 10 and 8 inches long. Take the longer piece and weave it through the aligned slot in the base as shown in the picture, with the soft side facing away from the base.



Wearable Clap Switch

ASSEMBLY GUIDE

Step 6


Take the smaller piece of hook and loop fastener and weave it through the angled slot with the soft side facing away from the base. Align it as shown in the picture.



Wearable Clap Switch ASSEMBLY GUIDE

Step 7

Place the battery in the base and route the wires as shown in the picture. Secure the switch in the base with a dab of CA glue.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 8

Connect the Circuit Playground to the computer and upload the `central_bleuartInstrucable0424.ino` file to it.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 9

Place the circuit playground on the four posts around the battery holder and use M3 screws to attach it onto each post.



Wearable Clap Switch ASSEMBLY GUIDE

Step 10

Take the lid and place it upside down. Insert both covers into the slots by the ports, and both buttons in the holes in the centre.



Wearable Clap Switch ASSEMBLY GUIDE

Step 11

Place the base assembly onto the lid, using the aligning peg to get the orientation correct. Secure the base onto the lid with two M3 screws.




Wearable Clap Switch

ASSEMBLY GUIDE

Receiver



Step 1

Cut two 6cm pieces of wire and strip 1cm from each end. Solder one end of each wire to the 3.5mm mono jack as shown in the picture. Connect the other ends to the NO and COM terminals on the relay.



Step 2

Cut 2 lengths of female headers to match the pins on the ItsyBitsy microcontroller. Attach the headers to the ItsyBitsy and place the headers onto the protoboard as shown in the picture. Solder the headers into place.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 3

Cut a 3 pin long section of female header and solder it to the protoboard as shown in the photo. Insert the battery backpack.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 4

Cut 4 10cm pieces of wire and strip 1cm from each end. Use two of these pieces to connect the 5V pin on the battery backpack header to the USB pin on the ItsyBitsy header, and the G pin on the battery backpack header to the G pin on the ItsyBitsy header



Wearable Clap Switch

ASSEMBLY GUIDE

Step 5

Use the remaining two wires and solder them to the centre and right pin of the toggle switch. The right pin is the one alignment pin is on, see this picture [more here]. Solder the other two ends to connect the BAT pin on the battery pack header to the BAT pin on the ItsyBitsy header.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 6

Cut 6 15cm pieces of wire and solder the end of each to every pin on the two potentiometers as in the attached photo. Take the right wires on both potentiometers and connect them to the G pin on the ItsyBitsy header, and both left wires to the 3V pin on the ItsyBitsy header. Connect the middle wire on one potentiometer to the A1 pin on the ItsyBitsy header, and the other to the A2 pin.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 7

Using breadboard wire, connect the – pin of the relay to the G pin of the ItsyBitsy header, the + pin to the battery pin, and the S pin to the 9 pin.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 8

Plug the battery into the battery backpack and screw the protoboard into the base using 4 M2 screws in the orientation shown in the picture.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 9

Remove the retaining nuts on both potentiometers and the toggle switch. Insert the potentiometer connected to A2 through the hole labeled Duration, orienting it with the tab, and securing it in place with the retaining nut. Repeat this with the potentiometer connected to A1 and the Delay hole, and the toggle switch and the ON/OFF hole.




Wearable Clap Switch

ASSEMBLY GUIDE

Step 10

Remove the retaining nut from the 3.5mm jack. Insert it through the last hole in the base, aligning the relay with the four posts on the base. Use two M3 screws on the rear posts to hold it in place. Secure the jack in place with the retaining nut.



Wearable Clap Switch

ASSEMBLY GUIDE

Step 11

Tuck the battery into the case and secure the lid with 4 M2 screws.

